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Ersatz-Behälter/ spare pot
klar/ transparent 1 St./ pc.

REF 07196KLAR

Pressure pot is independent of compressed air. No electrical parts. The working pressure of max. 2 bar is generated manually by turning the safety screw downwards giving pressure to the lid and therefore to the water.

- For the polymerisation of all self-curing resins up to 50 °C (122 °F) (e.g. direct relinings)
- Made of compact, thermostable plastic

TECHNICAL DATA

internal diameter of pot:	approx. 100 mm
height of useful space:	approx. 65 mm
gross weight:	approx. 1.3 kg

According to experiences of polymerisation	
water temperature:	50 °C (122 °F)
pressure load:	2 bar (29 PSI)
Polymerization time:	12-15 min.

HANDLING INSTRUCTIONS

1. Let the model thoroughly soak in water to prevent air bubbles.
2. Apply the resin and place the work into the pot.
3. Make sure that the air removal valve is closed.
4. Fill the pot with adequate hot water until overflowing and close the pot with the lid. A displacer at the bottom of the lid will push away superfluous water. This overflow avoids air being enclosed into the pot.
5. Insert the closed pressure pot into the pressure frame right back to the stop. Apply the necessary pressure by tightening the screw (maximum pressure load 2 bar).
6. After curing (time depends on the resin used) open the pressure screw and take the pot out of the pressure frame. To open the pot, equalize the pressure by turning the release valve (1/2 turn should be sufficient). Now we recommend to close the air removal valve.

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KENTZLER-KASCHNER DENTAL GmbH